

Ingenzi Tech | Admin Guide

This document serves as a comprehensive guide for administrators managing the LLM-Powered Philips HDI 5000 Ultrasound Machine Troubleshooting Tool. Follow the steps and guidelines carefully to configure, maintain, and troubleshoot the system effectively.

Table of Contents

1. **System Overview**
 2. **Prerequisites**
 3. **Setup and Configuration**
 - Environment Setup
 - API Key Management
 - Server Configuration
 4. **Running the Application**
 - React Frontend
 - Flask Backend
 - Chainlit Copilot
 - Feedback App
 5. **Maintenance and Monitoring**
 6. **Troubleshooting Common Issues**
 7. **Best Practices**
-

1. System Overview

The system is a multi-component application designed to assist with troubleshooting the Philips HDI 5000 ultrasound machine. The components include:

- **React Frontend:** A dashboard for users to interact with the system.
 - **Flask Backend:** Handles API requests and backend logic.
 - **Chainlit Copilot:** An AI-powered assistant accessible via a widget or standalone app.
 - **Feedback App:** Allows admins to manage and analyze user feedback.
-

2. Prerequisites

Ensure the following are available and installed:

1. **Software Requirements:**
 - Python (v3.8 or higher)
 - Node.js and npm
 - Chainlit (latest version)
 - Flask and required dependencies (via `requirements.txt`)
 2. **Hardware Requirements:**
 - Minimum 8GB RAM
 - At least 2 CPU cores
 - Stable internet connection for API integrations
 3. **Accounts and API Access:**
 - Google Cloud account with Calendar API enabled
 - GitHub account for managing repositories
 - OpenAI and Tavily accounts for API keys
-

3. Setup and Configuration

Environment Setup

1. Clone the repository:

```
Unset
git clone https://github.com/ahmedinhotahtiru/ingenzi_tech.git
cd ingenzi_tech
```

2. Create and activate python virtual environment

```
Unset
python -m venv capstone_venv
capstone_venv\Scripts\Activate.ps1
```

3. Install Python dependencies:

Unset

```
pip install -r requirements.txt
```

4. Install Node.js dependencies for the frontend:

Unset

```
cd ultrasound-dashboard
```

```
npm install
```

API Key Management

The system requires various API keys for external services. Follow these steps:

1. **Google Calendar API:**

- Navigate to the [Google Cloud Console](#).
- Enable the Calendar API and generate OAuth 2.0 Client IDs.
- Replace placeholders in `ultrasound-dashboard/Home.js`:

Unset

```
const CLIENT_ID = 'your-client-id.apps.googleusercontent.com';  
  
const API_KEY = 'your-api-key';
```

2. **GitHub Token:**

- Create a personal access token from [GitHub Settings](#).
- Replace the token in `backend/end_points.py`:

Unset

```
g = Github("your-personal-access-token")
```

3. OpenAI and Tavily API Keys:

- Generate keys from [OpenAI](#) and [Tavily](#).
- Replace placeholders in `chat_ultrasound_chroma.py`:

Unset

```
os.environ["TAVILY_API_KEY"] = 'your-tavily-api-key'  
  
embeddings =  
OpenAIEmbeddings(openai_api_key="your-openai-api-key")
```

Server Configuration

If any port is occupied, manually update the port configurations:

- React Frontend: Modify the `npm start` command.
- Flask Backend: Update the `app.run` port in `backend/end_points.py`.
- Chainlit Copilot: Adjust the `chainlit run` command.

4. Running the Application

1. Start the React Frontend:

- Navigate to `ultrasound-dashboard` and run:

Unset

```
npm start
```

Screenshot of running app below:

```
Windows PowerShell
Compiled successfully!

You can now view ultrasound-dashboard in the browser.

Local:      http://localhost:3000
On Your Network: http://192.168.1.70:3000

Note that the development build is not optimized.
To create a production build, use npm run build.

webpack compiled successfully
```

2. Start the Flask Backend:

- Navigate to `backend` and run:

Unset

```
python .\end_points.py
```

```
Windows PowerShell
(capstone_env) (base) PS C:\Users\Ahmed Issah Tahiru\Desktop\CMU-Africa\Second Year\EAI Project Methods\Copilot LLM\back
end> python .\end_points.py
* Serving Flask app 'end_points'
* Debug mode: on
WARNING: This is a development server. Do not use it in a production deployment. Use a production WSGI server instead.
* Running on http://127.0.0.1:5000
Press CTRL+C to quit
* Restarting with stat
* Debugger is active!
* Debugger PIN: 941-580-033
127.0.0.1 - - [13/Dec/2024 23:39:51] "GET /api/last-service-date HTTP/1.1" 200 -
127.0.0.1 - - [13/Dec/2024 23:39:51] "GET /api/last-service-date HTTP/1.1" 200 -
```

3. Start the Chainlit Copilot:

- Navigate to **RAG** and run:


Unset

```
chainlit run .\chat_ultrasound_chroma.py
```

```
Windows PowerShell
(capstone_env) (base) PS C:\Users\Ahmed Issah Tahiru\Desktop\CMU-Africa\Second Year\EAI Project Methods\Copilot LLM\RAG>
chainlit run .\chat_ultrasound_chroma.py
2024-12-13 23:40:41 - Anonymized telemetry enabled. See https://docs.trychroma.com/telemetry for more information.
C:\Users\Ahmed Issah Tahiru\Desktop\CMU-Africa\Second Year\EAI Project Methods\Copilot LLM\RAG>.\chat_ultrasound_chroma.py:164: LangChainDeprecationWarning: The class 'Chroma' was deprecated in LangChain 0.2.9 and will be removed in 1.0. An updated version of the class exists in the :class:'~langchain-chroma' package and should be used instead. To use it run 'pip install -U :class:'~langchain-chroma' and import as 'from :class:'~langchain-chroma import Chroma''.
  vector_store = Chroma(

Loaded existing collection: ultrasound_manuals
Loaded existing collection: error_manuals
Loaded existing collection: maintenance_manuals
Loaded existing collection: device_history
2024-12-13 23:40:49 - Your app is available at http://localhost:8000
2024-12-13 23:40:55 - Translated markdown file for en-US not found. Defaulting to chainlit.md.
```

The copilot now runs on top of the react application as a chat widget shown below:



4. Start the Feedback App:

- Navigate to **RAG** and run:

Unset




```
python feedback_app.py
```

Screenshot of running feedback shown below:

```
Windows PowerShell
[capstone_env] (base) PS C:\Users\Ahmed Issah Tahiru\Desktop\CMU-Africa\Second Year\EAI Project Methods\Copilot LLM\RAG>
python .\feedbackash.py
Dash is running on http://127.0.0.1:8050/

* Serving Flask app 'feedbackash'
* Debug mode: off
WARNING: This is a development server. Do not use it in a production deployment. Use a production WSGI server instead.
* Running on http://127.0.0.1:8050
Press CTRL+C to quit
```

The running feedback app looks like this:

Ingenzi Feedback			
ID	User Prompt	Feedback	Helpful
11c1bbc3-c77a-4dfc-8bb9-b67c1ae561e3	i have error 0060	This is the correct error code description for error 0060	
2b7599ef-409c-4b24-af3d-db8dbf93c0ba	the image is not clear when i scan	This was a great resolution to improving scan quality	
cc954d17-5243-4b29-9c23-2d1b650056ff	retrieve device logs	Worked great	
7f6c3986-d161-4ce1-839a-22f19902272b	give me best practices	Too generic response	

5. Maintenance and Monitoring

1. **Logs:**
 - Monitor logs for each component to identify issues.
 - Flask logs: `backend/logs`
 - React logs: Displayed in the terminal during `npm start`.
2. **API Quotas:**
 - Monitor usage limits for Google, OpenAI, and Tavily APIs.
3. **Updating Dependencies:**
 - Periodically update Python and Node.js dependencies:

Unset

```
pip install --upgrade -r requirements.txt
```

```
npm update
```

6. Troubleshooting Common Issues

Issue	Possible Cause	Solution
Frontend not loading	React server not running	Ensure <code>npm start</code> is executed.
Backend API not responding	Flask server not running	Check logs and rerun <code>end_points.py</code> .
Chainlit copilot not accessible	Port conflict	Change the port and rerun.
Feedback app not opening	Port <code>8050</code> is occupied	Change the port in <code>feedback_app.py</code> .
API key authorization errors	Invalid or expired keys	Regenerate keys and update configuration files.

7. Best Practices

1. Security:

- Store sensitive keys in environment variables or `.env` files.
- Avoid committing API keys to version control.

2. Backups:

- Regularly back up critical files and databases.

3. Documentation:

- Maintain updated documentation for all custom configurations.